drive

DeVilbiss
HEALTHCARE

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Instructions for use

CONTENTS

Preface	3
Introduction.....	3
Declaration of conformity	3
Indication/Contraindication.....	3
Scope of delivery.....	3
Safety instructions.....	4
Electromagnetic compatibility (EMC)	4
Setting up the system.....	4
Explanation of functions.....	5
Startup	6
Hand check.....	6
Alternating pressure function.....	6
Static function.....	6
Autofirm function	6
Seat function.....	6
Air flow	7
CPR (cardiopulmonary resuscitation).....	7
Transport function	7
Power failure	7
Alarm function.....	7
Cleaning, disinfection and care	8
Reuse/patient change	8
Maintenance, service and storage.....	9
Disposal	9
Warranty.....	9
Troubleshooting.....	10
Specifications	10
Type label.....	11
Service life	11

Preface

Thank you for choosing an alternating pressure system from DRIVE MEDICAL.

Please read these instructions for use carefully prior to initial use. They contain important safety information and valuable tips for the proper use and care of the Med Aire FlexWave alternating pressure system. If you have questions or need additional information, contact your DRIVE MEDICAL specialist dealer who supplied your alternating pressure system.

All repairs as well as some adjustments require special technical training and must therefore be made by the medical supplies dealer.

Introduction

The Drive Medical alternating pressure system Med Aire FlexWave is specially developed for the prophylaxis and treatment of decubitus, up to and including grade 4 decubitus according to EPUAP. This system replaces the existing mattress in your bed.

Equipped with a cell-on-cell system and a high-performance unit with an automatic start-stop mechanism, this system is ideally suited for long-term therapy and facilitates home care.

The Med Aire FlexWave can be operated in alternating pressure mode as well as in static mode. Through an integrated seat function, the pressure in the cells can be adapted to a seated position. A special autofirm function facilitates the mobilization of the patient on this system.

The system is suitable for patients from 20 kg to a max. of 200 kg.

IMPORTANT! The Med Aire FlexWave does not replace regular repositioning of the patient!

Declaration of conformity

For the alternating pressure system Med Aire FlexWave, DRIVE MEDICAL GMBH & Co. KG declares compliance according to the EC directive 93/42 for medical devices, taking amending directive 2007/47/EC into account.

Indication/Contraindication

Indication:

- Decubitus prophylaxis for patients at high to very high risk
- Pressure distribution and relief for immobile and partially mobile patients
- Decubitus therapy up to and including grade 4 according to EPUAP

Static pressure-relief positioning is suitable for:

- Patients who only occasionally tolerate dynamic alternating pressure or who do not tolerate it at all.
- Pain patients

Contraindication:

- Body weight less than 20 kg and over 200 kg
- In the case of unstable fractures (particularly of the back, cervical and lumbar spine)
- In the case of neurological diseases in which dynamic alternating pressure should not be used, only static pressure-relief positioning should be used.

To avoid uncertainties regarding the indication and contraindication, we recommend consulting the attending physician.

Scope of delivery

Please check the contents for damage and completeness. If one of the parts listed below is missing, please contact your specialist dealer.

1. Med Aire FlexWave unit
2. Med Aire FlexWave mattress replacement system (including cell-on-cell system, connection hose)
3. Breathable, water-repellent and quilted cover
4. Instructions for use

Safety instructions

- Keep the unit away from water and/or other liquids.
- If water or another liquid gets into the device, disconnect it from the mains immediately.
- Do not allow the system to come into contact with open flames or smoldering objects.
- Do not operate the unit in the vicinity of flammable substances or chemicals => Risk of explosion!!!
- The pump housing may only be opened by authorized technical personnel => Risk of electric shock!!!
- Avoid heavy strain and/or damage to the housing and/or power cord.
- Avoid pulling, kinking and running over the power cord.
- If the device fuses are defective, please have the device checked by an authorized specialist dealer.
- Use this product only for the purpose described in the instructions for use.
- Operate the device only with the specified voltage (230 V/50 Hz).
- The unit should not be covered => Risk of overheating!
- Connect the device only to the mattress intended for it.
- This device is safeguarded against electrical shock according to the specifications for type BF devices! Product protection class I
- Do not allow the system to come into contact with sharp objects (knives, scissors, cannulas, etc.).
- This device should not be disposed of with household waste. More information is available from local waste management companies.
- **When using side rails on adjustable beds for disabled persons, a sufficient distance between the top edge of the positioning system and the upper end of the bed rail must be ensured. (220 mm according to DIN 1970:2000). If necessary, corresponding side rail risers should/must be attached to protect the patient.**

Electromagnetic compatibility (EMC)

This product meets the protection requirements listed in the council directive on the approximation of the laws of the member states on electromagnetic compatibility (EMC, IEC 60601-1-2). However, mutual interference between electrical devices, particularly when mobile telephones are used, cannot be completely ruled out under certain conditions.

Setting up the system

The Med Aire FlexWave replaces the existing mattress in your standard or hospital bed.

It does not impair the functionality of the standard or hospital bed.

- Take the Med Aire FlexWave out of the packaging.
- Check the contents for completeness and any visible damage.
- Now remove the existing mattress from the standard or hospital bed if you have not already done so.
- Place the mattress replacement system in the existing standard or hospital bed. The foot end is marked with 2 feet.
- Secure the mattress replacement system using the available straps to the movable parts of the slatted frame to avoid impairing the functionality of the system and the adjustable slatted frame.
- Place the unit using the hooks provided at the foot end, on the floor, or on an even surface.

Do not in any way cover the unit. Risk of overheating!!!

Now attach the connection hose to the unit.

Important! The connection hose must audibly click into place in order to prevent an unintentional separation between the unit and the mattress! Kinking or twisting the connection hose must be avoided.

Now the Med Aire FlexWave is ready for use.

Explanation of functions

Weight setting

User weight setting

Autofirm

Maximum increase of cell pressure in all cells

Static and alternating pressure

Adjustment of static pressure-relief positioning and dynamic alternating pressure.

Cycle time

Adjustable: 10/15/20/25 minutes

Lock

Automatically enabled after 30 seconds of inactivity or by pressing the "Tastensperre/Lock" button.

Mute

Acoustic alarm signal can be muted and remains muted until the button is pressed once again!

Seat

Increase in the cell pressure in a seated position

Startup

Connect the plug to the unit and an appropriate power source. (See Fig. 1).

Switch on the unit. The on and off switch is located on the side of the device, above the power cord. (See Fig. 2).

The mattress is now in the filling phase while the LED "Low Pressure" is lit. (See Fig. 3).


Adjust the weight of the patient using the (+) and (-) button on the weight adjustment scale*; if the "Low Pressure" LED goes out and the "Normal Pressure" LED is lit, the patient can be placed on the system.

***The weight adjustment scale indicated is a guideline and does not reflect the actual pressure in the air chambers.**

Hand check

After each change in pressure and position, use your hand to check the proper pressure setting at the buttocks/hips. When correctly adjusted, a hand should be able to be easily inserted between the patient and a deflated cell.

Alternating pressure function

All air chambers are evenly inflated and deflated in a 1:1 ratio. An alternating cycle involves in each case the inflating and deflating of an air chamber and takes approx. 10/15/20/25 minutes, depending on the setting. The first 3 air chambers in the head area remain static.

The "Alternating Pressure" LED indicates that the alternating pressure function is enabled.

Static function

The Med Aire FlexWave has a static function as a standard feature.

This is switched on using the button. The "Static" LED indicates that the static function is enabled.

The air chambers all maintain the same pressure and ensure adjustable pressure-relief positioning.

In static pressure-relief mode, the deepest point on the patient (buttock/hip area) should sink into the system to approximately halfway down the upper cell system (approx. 3-4 cm).

Note: Perform a hand check to verify; see hand check section.

Autofirm function

This function facilitates the positioning and mobilization of the patient through a firm and stable surface. The pressure within the cells is increased to a maximum level. For reasons of safety, this function automatically switches off after 30 minutes!

Seat function

In this case, the pressure in the air chambers is increased and adapted to the short-term seated position. The objective is to prevent the patient from "bottoming out."

The seat function can be manually switched on and off. **This function may only be used for short periods of time.**

For safety reasons, this function automatically switches off after 90 minutes.

Important!: The system is not suitable for long periods of sitting in bed!

Air flow

The Med Aire FlexWave has nine air flow cells. This group of air chambers is microperforated at the edge and in this way, it ensures that the Med Aire FlexWave is actively supplied with air.

CPR (cardiopulmonary resuscitation)

Before resuscitation can/may be performed on the system, the system must be fully deflated.

For rapid deflation, pull down on the red tab printed with "CPR" at the head end (see fig.). The tab is on the same side as the hose system. The valves open and air escapes from all of the cells.

To supply the system with air once again, close the CPR valve again by pressing it firmly and securely fasten the tab marked "CPR" once again.

Transport function

Remove the hose system from the unit and promptly close it with the attached cap. Now the mattress is in transport mode. To return to normal operation, remove the cap on the connection hose and connect it back to the unit and start the unit.

Note: When the Med Aire FlexWave is in transport mode, it is possible to leave patients on the system for up to 12 hours. However, at this time, we recommend regularly monitoring the patient and the mattress replacement system.

Power failure

In the event of a power failure, set the system to transport mode. See section on transport function.

Alarm function

Low pressure

The "Low Pressure" LED is lit whenever the unit detects too little air in the cells. An acoustic signal sounds at the same time if it has not been switched off beforehand by pressing the "Mute" button.

Once the set pressure has been reached in the cells again, the display goes out and the acoustic signal stops. The "Normal Pressure" LED is lit.

Important: During the filling phase, the "Low Pressure" LED is lit; once the pressure in the air chambers has reached the set user weight, the "Low Pressure" LED goes out.

System failure/service display

If the unit detects an internal error, the "System Failure" LED is continuously lit and an acoustic signal additionally sounds. This can be muted using the "Mute" button. The "System Failure" LED will still be lit. Please notify your specialist dealer.

Power failure

If there is a power failure or if the unit is disconnected from the power supply, the "Power Failure" LED begins to blink and an acoustic signal sounds.

This alarm stays on until the unit is reconnected to the mains or the On/Off switch is operated.

Cleaning, disinfection and care

Unit

Before cleaning, disconnect the unit from the power supply. During cleaning, please ensure that no fluids/cleaning solutions come into contact with the plug.

The unit should not be saturated or soaked.

The unit may be cleaned by hand. To do this, use a moist (not wet) cloth and a mild household cleaning and disinfectant solution.

Do not use any solutions which contain phenol or abrasive cleaners. This can corrode the surface of the unit.

Allow the unit to dry thoroughly.

Perform a function test after the device is thoroughly dry.

Mattress replacement system

The mattress replacement system should not be saturated or soaked*, however it can be cleaned by hand. To do this, use a moist (not wet) cloth and a mild standard cleaning and disinfectant solution.

*** There is a risk of liquid getting into the cell chambers and hoses.**

Cover

The cover can be cleaned by hand. To do this, use a moist (not wet) cloth and a mild cleaning solution.

You can detach the cover using the zipper and wash it at a max. temperature of 60° using a household detergent.

Do not dry the mattress cover in the dryer.

Note! If the Med Aire FlexWave is machine-washed and -disinfected, a validated and automated method according to the Robert Koch Institute (RKI) and/or the Association for Applied Hygiene (VAH) is necessary.

Reuse/patient change

In the case of reuse/patient change, the Med Aire FlexWave must be disinfected and cleaned according to a validated method recognized by the Robert Koch Institute.

During processing, take note of the medical device law, the "Guidelines for the Processing of Medical Devices" as well as the "Guidelines for Hospital Hygiene and Infection Prevention" of the Robert Koch Institute.

Maintenance, service and storage

Maintenance

This product is subject to the Medical Device Law (MPG) and the Medical Devices Operator Ordinance (MPBetreibV). Internationally, the respective valid national regulations apply.

The functionality and proper condition must be checked before each use.

The following points must be regularly checked:

- Check power cord and plug for wear and damage.
- Check the housing of the unit for damage.
- Check the alternation between air supply and ventilation in each air chamber circuit. This alternation must take place according to the cycle time of 10/25/20/25 minutes set.
- Check all hoses and connections for leaks, kinks and breakage.
- Check all LEDs which must light up in their respective function.

Please contact your specialist dealer if you identify any damage or defects.

During all visual inspections and maintenance, the unit must be disconnected from the power supply.

Testing and assessments must be performed and documented by suitable technical staff in accordance with DIN EN 62353.

Service

Only original replacement parts and accessories from Drive Medical may be used. If the system develops a defect or fault, please contact a trusted specialist dealer.

Med Aire FlexWave can be checked and repaired by an authorized specialist dealer or via your specialist dealer from Drive Medical.

Drive Medical recommends performing a technical safety check at least every 2 years.

Storage

- Remove the hose system from the compressor.
- Fully deflate the mattress.
- After you have fully deflated the overlay, roll up the overlay, starting at the foot end. You can use the fastening strap to secure the rolled-up overlay to prevent the Med Aire FlexWave from unrolling.

Disposal

Please contact your local waste management company regarding professional disposal.

Please be aware that the alternating pressure overlay and the mattress cover may be carrying microbes. This can lead to transmission of infection. Please dispose of the overlay and mattress cover such that there is no risk to you and third parties.

Warranty

Drive Medical grants a 2-year warranty for the product against material and manufacturing errors.

Damage due to natural wear and/or improper handling remain excluded from the warranty.

Damage due to material and/or manufacturing errors will be corrected free of charge through a replacement or repair.


Troubleshooting

Symptom	Remedy
Unit does not work	<ul style="list-style-type: none"> Check whether the plug is correctly inserted. Switch on the unit again. If the LED on the On/Off switch does not light up, check the voltage of the outlet or insert the plug into another outlet. If the LED on the On/Off switch lights up but the unit does not work, please contact your specialist dealer.
"Low Pressure" alarm	<ul style="list-style-type: none"> Check whether the connection hose on the unit is correctly inserted. The CPR valve must be closed.
Mattress does not pump up	<ul style="list-style-type: none"> Ensure that all cells are correctly connected to each other and that the connection hoses are not twisted or kinked.
Patient bottoms out	<ul style="list-style-type: none"> This indicates excessively low pressure in the air chambers. Increase the user weight and perform a hand check to check whether the patient is sinking down. Check air filter (located above the on/off switch on the side of the unit); please clean or replace it. If the unit continues to generate too little pressure, the unit must be replaced. Contact your specialist dealer
No alternating pressure takes place	<ul style="list-style-type: none"> Ensure that the air supply hoses are not kinked or twisted. The static function must be switched off. When the static function is switched on, no alternating pressure takes place. Ensure that the autofirm and seat functions are switched off.


Specifications

Unit	Mattress replacement system		Environmental parameters
Power supply	230 V / 50 Hz Protection class 1; IPX0; Type BF, AP/APG NO	Dimensions 200 x 88 x 20 cm	Temperatures Operation 10 – 35°C Storage -15 – 50°C Shipment -15 – 70°C
Power consumption	Max. 25 watts	Cells Nylon/PU; 20 cells of which 9 are air flow cells, 3 are static head cells	Relative humidity Operation 20 – 80% => non-condensing Storage 10 – 90% => non-condensing
Weight	3.9 kg	Cell height 20 cm	
Dimensions	36.9 x 13.5 x 25 cm	Cover Nylon/PU, breathable, washable	
Pressure range	10- 50 mbar		
Cycle time	10/15/20/25 min.		
Max. load	20 - 200 kg		

Type label

1. Product name
2. Article number
3. Serial number
4. Production date year/month
5. Do not dispose of in household waste
6. Follow the instructions for use
7. CE marking
8. Do not dry in dryer
9. Do not iron
10. Wash in a normal wash cycle up to a temperature of 60°C
11. Do not bleach

1. Product name
2. Article number
3. Serial number
4. Production date year/month
5. Voltage
6. DIN standard
7. Follow the instructions for use
8. CE marking
9. Type BF
10. Protection type
11. Protection class 1 and 2
12. Do not dispose of in household waste
13. Important

Service life

Our company assumes a long service life for this product as long as the product is used as intended and all maintenance and service requirements are observed.

This service life can be significantly prolonged if the product is carefully handled, maintained, cared for and used.

The service life may also be significantly shortened through extreme and improper use, however.

Determination of the service life by our company does not represent any additional guarantee.

Your specialist dealer:

Subject to error and change
Last updated: 12.21.2017

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